

		Inner Diameter		Guidewire compatibility	Microcatheters													
					DraKon™						SeQure®							
					1.7Fr	1.9Fr	2.4Fr	2.7Fr	2.8Fr	3.0Fr	1.7Fr	1.9Fr	2.4Fr	2.7Fr	2.8Fr	3.0Fr		
Angiographic Catheters	4Fr	Terumo Glidecath®	0.041"	1.033mm	0.038"	✓	✓	✓	✓			✓	✓	✓	✓			
		Cordis Tempo®	0.042"	1.067mm	0.038"	✓	✓	✓	✓			✓	✓	✓	✓			
		Merit Impress®	0.040"	1.02mm	0.038"	✓	✓					✓	✓					
		Merit Impress®	0.040"	1.02mm	0.035"													
		BSC Imager™ II Selective	0.039"	0.99mm	0.035"													
	5Fr	Terumo Glidecath®	0.043"	1.10mm	0.038"	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Terumo Glidecath® XP	0.047"	1.20mm	0.038"	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
		Cordis Tempo®	0.043"	1.09mm	0.038"	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓		
		Merit Impress®	0.046"	1.07mm	0.038"	✓	✓	✓	✓			✓	✓	✓	✓			
		Merit Impress®	0.046"	1.07mm	0.035"													
BSC Imager™ II Selective		0.045"	1.14mm	0.035"														
BSC Imager™ II Selective		0.050"	1.27mm	0.038"	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

DraKon™ and SeQure® microcatheters are class IIb medical devices intended for use by interventional radiologists and interventional oncologists for the infusion of contrast media into all peripheral vessels and for drug infusion in intra-arterial therapy, and infusion of embolic materials. They should not be used in cerebral vessels.

For complete information about precautions and optimal usage conditions for these medical devices, we recommend consulting the instructions for use supplied with each device or by your local Guerbet representative(s). This information is intended for use only in countries with appropriate regulatory authority market clearance. Notified Body: MedCert 0482.

Manufacturer: Accurate Medical Therapeutics Ltd. EC Rep: Guerbet. Document creation date: December 2020.

DraKon™ and SeQure® are a registered trademark of Guerbet Group or its affiliates.

Glidecath® and Glidecath® XP are a registered trademark of Terumo Group.

Tempo® is a registered trademark of Cordis Group.

Impress® is a registered trademark of Merit Group.

Imager™ II Select is a registered trademark of Boston Scientific Group.

P20000407